

THEME: INTEGRITY AND SECURITY IN THE MANAGEMENT OF EXAMINATION IS MY RESPONSIBILITY

Candidate's Name

NAHWERA BESUMA SETH

Signature

*[Handwritten Signature]*

Random No.

Personal No.

0 8 2

(Do not write your School/Centre Name or Number anywhere on this booklet)

545/3

CHEMISTRY  
(PRACTICAL)

Oct./Nov. 2019

2 hours



UGANDA NATIONAL EXAMINATIONS BOARD

Uganda Certificate of Education

CHEMISTRY PRACTICAL

Paper 3

2 hours

### INSTRUCTIONS TO CANDIDATES:

Answer **both** questions. Answers are to be written in the spaces provided in this booklet.

All your work **must** be in **blue or black ink**. Any work done in pencil **except** drawings will **not** be marked.

You are **not** allowed to use any reference books (i.e. text books, booklets on qualitative analysis etc.).

All working must be clearly shown.

Mathematical tables and silent non-programmable calculators may be used.

For Examiners' Use Only		
Q.1	25	
Q.2	25	
Total	50	

1. You are provided with the following:

BA1, which is a solution made by dissolving 3.45 g of a hydrated salt  $X.nH_2O$  in  $250\text{ cm}^3$  of water.

BA2, which is a 0.1 M hydrochloric acid.

You are required to determine the value of  $n$  in the salt.

### Procedure:

Pipette  $25\text{ cm}^3$  (or  $20\text{ cm}^3$ ) of BA1 into a conical flask. Add 2-3 drops of methyl orange indicator and titrate with BA2 from the burette.

Repeat the titration until you obtain consistent results.

Record your results in the table below.

Volume of pipette used...  $25.0$   $\text{cm}^3$   $\frac{1}{2}$  ( $\frac{1}{2}$  mark)

0. decimal  
one decimal

Final burette reading  
to correct if  
it is in range  
of Centre average

Final burette reading ( $\text{cm}^3$ )	29.10	28.00	29.60
Initial burette reading ( $\text{cm}^3$ )	1.50	0.50	2.10
Volume of BA2 used ( $\text{cm}^3$ )	27.60	27.50	27.50

7 1/2

Final burette reading to two decimal places 1 mark @  
Initial burette reading to two decimal places 1 mark @  
(7 1/2 marks)

Titre values of BA2 used for average volume of BA2 used  $\frac{1}{2}$  (01 mark)

27.50 and 27.50  $\text{cm}^3$  1 mark

Average volume of BA2 used

$$\frac{27.50 + 27.50}{2} = 27.50 \pm 0.1 \pm 0.2 \pm 0.3 \pm 0.4 \pm 0.5 \text{ cm}^3$$

(2 1/2 marks)

### Questions

(a) Calculate the;

(i) number of moles of hydrochloric acid that reacted. ( $1\frac{1}{2}$  marks)

1000  $\text{cm}^3$  of hydrochloric Acid contain 0.1 Mole

27.50  $\text{cm}^3$  of hydrochloric Acid contain  $0.1 \times \frac{27.50}{1000}$  moles  
 $= 0.00275$  moles  
At least 4 d. places



- (ii) number of moles of  $X.nH_2O$  that reacted. (1 mole of  $X.nH_2O$  reacts with 2 moles of hydrochloric acid). (02 marks)

1 mole of  $X.nH_2O$  reacts with 2 moles of hydrochloric acid

Moles of  $X.nH_2O$  that reacted  $\frac{1}{2} \times 0.00275$  ✓  
 $= 0.001375$  moles ✓ 02

- (iii) number of moles of  $X.nH_2O$  in  $250 \text{ cm}^3$  of BA1. (03 marks)

$25 \text{ cm}^3$  of BA<sub>1</sub> contain 0.001375 moles of  $X.nH_2O$

$\therefore 250 \text{ cm}^3$  of BA<sub>1</sub> contain  $\frac{0.001375 \times 250}{25}$  moles ✓ 03  
 $= 0.01375$  moles

- (b) Determine the value of  $n$  in  $X.nH_2O$ . (5½ marks)  
 (H = 1; O = 16; X = 106)

0.01375 moles of  $X.nH_2O$  weigh 3.45g ✓

$\therefore$  1 mole of  $X.nH_2O$  weighs  $\frac{3.45 \times 1000}{0.01375} = 250.91 \text{ g}$  ✓

Mass of  $nH_2O = 250.91 - 106 = 144.91$  ✓

Rmm of  $H_2O = (1 \times 2) + 16 = 18$  ✓

$\therefore 18n = 144.91$

$n = \frac{144.91}{18} = 8$  05½

$n = 8$  ✓

25

Alternative method

moles in  $1000 \text{ cm}^3$

moles in  $250 \text{ cm}^3$  (missing part)



2. You are provided with substance Q, which contains two cations and one anion. Carry out the following tests to identify the cations and the anion present in Q. Identify any gas(es) that may be evolved. Record your observations and deductions in the table below.

(25 marks)

TESTS	OBSERVATIONS	DEDUCTIONS
(a) Dissolve one spatula end-ful of Q in about 5 cm <sup>3</sup> of water. Add excess sodium hydroxide solution; shake well and filter. Keep both the filtrate and the residues.	Brown solution ✓ Brown Residue ✓ or brown ppt insoluble Colourless filtrate ✓	Fe <sup>3+</sup> ✓ Al <sup>3+</sup> ✓ Pb <sup>2+</sup> ✓ Zn <sup>2+</sup> ✓ Ignore NH <sub>4</sub> <sup>+</sup>
(b) To the filtrate, add dilute nitric acid until the solution is just acidic. Divide the acidified solution into five portions.	White ppt soluble in excess ✓	Al <sup>3+</sup> ✓ Pb <sup>2+</sup> ✓ Zn <sup>2+</sup> ✓ Al <sup>3+</sup> ✓ Pb <sup>2+</sup> ✓ Zn <sup>2+</sup> (Ignore)
(i) To the first portion of the acidified solution, add dilute sodium hydroxide drop-wise until in excess.	White ppt soluble in excess ✓	Al <sup>3+</sup> ✓ Pb <sup>2+</sup> ✓ Zn <sup>2+</sup> ✓
(ii) To the second portion of the acidified solution, add aqueous ammonia drop-wise until in excess.	White ppt insoluble in excess ✓	Al <sup>3+</sup> ✓ Pb <sup>2+</sup> ✓ but appeared in
(iii) To the third portion of the acidified solution, add potassium iodide solution.	No observable change ✓ No yellow ppt. Solution remained colourless ✓	Pb <sup>2+</sup> absent ✓ Al <sup>3+</sup> present ✓ as appeared in (i) above.

TESTS	OBSERVATIONS	DEDUCTIONS
(iv) To the <b>fourth</b> portion of the acidified solution, add lead(II) nitrate solution and warm.	White ppt Insoluble on warming	$\text{SO}_4^{2-}$ ✓ 0 1/2
(v) Use the <b>fifth</b> portion of the acidified solution to carry out a test of your own to confirm the anion in Q.  TEST Added $\text{Ba}(\text{NO}_3)_2(aq)$ ✓ OR $\text{BaCl}_2(aq)$	White ppt ✓	$\text{SO}_4^{2-}$ ✓ 03 2 1/2
(c) Dissolve the residue in minimum amount of dilute sulphuric acid and divide the resultant solution into two parts. (i) To the <b>first</b> part of the solution add sodium hydroxide solution drop-wise until in excess.	Pale Brown solution ✓ Acc. yellow soln ✓ Brown ppt Insoluble in excess ✓	$\text{Fe}^{3+}$ ✓ 02
(ii) To the <b>second</b> part of the solution, add 1 small piece of zinc granules and leave the solution to stand for 5 minutes. Divide the solution into two portions and use them for part (d).	Pale yellow solution turned green ✓ soln 2	$\text{Fe}^{3+}$ ✓ $\text{Fe}^{3+}$ ions reduced to $\text{Fe}^{2+}$ ions ✓ Accept $\text{Fe}^{2+}$ ions 0 1/2



TESTS	OBSERVATIONS	DEDUCTIONS
(d) (i) To the first part of the solution, add sodium hydroxide solution drop-wise until in excess.	Green ppt Insoluble in excess	$\text{Fe}^{2+}$ formed $0\frac{1}{2}$
(ii) To the second part of the solution, add aqueous ammonia drop-wise until in excess.	Green ppt Insoluble in excess	$\text{Fe}^{2+}$ formed $\checkmark$ $\therefore \text{Fe}^{3+}$ present in Q $0\frac{2}{2}$

- (c) (i) The cations in Q are  $\text{Al}^{3+}$  Correctly confirmed in b(iii)  $0\frac{1}{2}$
- (ii) The anion in Q is  $\text{Fe}^{3+}$  Correctly confirmed in c(i) and a, s(ii)  $0\frac{1}{2}$
- $\text{SO}_4^{2-}$  Correctly confirmed in b(v)  $0\frac{1}{2}$

Symbols MUST be correctly written.

An extra wrong ion cancels the correct one.

All technical words MUST be correctly written  
eg precipitate, Residue filtrate soluble, Insoluble

25

No final burette reading: Award for Initial  
Deny for Volume of B.A. used.  
Award for ~~Average~~ <sup>Range</sup> if ~~in the range~~  
Deny accuracy.

2 Columns filled: Mark normally

1 Column filled: Mark normally only the filled  
Column.

Initial greater than  
final reading: No mark.

No volume of pipette indicated

→ Deny for Final burette reading and Accuracy  
(Award the rest.)

→ Award full marks if volume of pipette is  
reflected in a(ii)